

## Autoclaving the InstaWarmer

The InstaWarmer pack can be reset and sterilized in a standard hospital autoclave cycle.

**Sterilization Time: 15 Minutes**

**Dry Time: 10 Minutes**

**Pressure: 100-105 Bar/PSI**

### Step 1:

Wipe InstaWarmer pack clean with alcohol wipes, or other sanitizing solution



### Step 2:

Roll up the pack before placing in autoclave. Please wrap in autoclave paper.



### Step 3:

Place the pack inside the autoclave. Close autoclave door and seal using standard procedures. Seal tight for good seal.



Step 4: Switch on the autoclave with gravity cycle at 105 degrees Celsius for 15 minutes total sterilization time. Use 100-105 bar/PSI. Close per standard procedure.



Step 5: Switch on the autoclave dry cycle at 105 degrees Celsius for 10 minutes.

